



GXR-S SERIES

High Performance and Reliability

GXR DIAGNOSTIC X-RAY SYSTEM

High Performance and Reliability

GXR Diagnostic X-ray System for Digital Era

The new designed GXR diagnostic X-ray system provides an analogue radiographic room that perfectly fits your workflow and budget, but it can be easily upgraded to DR system with the help of DR interface and PC interface in the generator as well as Bucky suitable to Flat Panel Detector.

DRGEM's GXR X-ray system is equipped with a high frequency X-ray generator which consistently produces high quality images based on high quality X-ray output with a very small kV ripple and accurate mA and timer control.

GXR X-ray system is designed with complete consideration for both operator and patient comfort, and specially designed for minimizing electrical lock noise.

Experience the quality, durable and user-friendly operation for all applications with the GXR X-ray system!



Floor-Ceiling Mounted Tube Stand




Four-way Tabletop Patient Table



Elevating Patient Table




Vertical Wall Bucky Stand




Floor Mounted Tube Stand

DR Upgrade Option

GXR System (Film)



DR System





© 2011 DRGEM Corporation - All rights reserved.

Reproduction in whole or in part is prohibited without the written consent of the copyright-holder.

Specifications and features are subject to change without notice.

All changes will be in compliance with regulations governing manufacture of medical equipment.

7th Fl. Dongil Techno-town 1st, 826 Guro-dong, Guro-gu, Seoul, Korea 152-879
Tel: +82-2-869-8566 / Fax: +82-2-869-8567 / www.drgem.co.kr / drgem@drgem.co.kr